


<b>Naziv:</b>	<b>Vickers microhardness device equipped with CCD Image Automatic Measuring System</b>
<b>Inventarni broj</b>	F 00517
<b>Tip opreme:</b>	Research and teaching equipment
<b>Lokacija:</b>	Physics department
<b>Dostupnost:</b>	Internally/Commercially
	
<b>Tehničke karakteristike:</b>	<p>Model: DHV-1000Z  Vickers microhardness testing  Load: 10 g, 25 g, 50 g, 100 g, 200 g, 300 g, 500 g, 1000 g  Objective: 10x i 40 x  Eyepiece: 10x  Loading time: 0-60 s  Maximum height of a sample: 75 mm  Maximum width of a sample: 95 mm  CCD camera: 1.30 million pixel</p>
<b>Godina proizvodnje:</b>	2020.
<b>Usluge:</b>	Microhardness measurements
<b>Radno vrijeme:</b>	8-16 h
<b>Kontakt:</b>	<p>Doc. dr. Amra Salčinović Fetić  033 279 953  amra.s@pmf.unsa.ba</p>